



BIO AND NEURO AMPLIFIERS

Owner's Guide

This document was, as far as possible, accurate at the time of release. However, changes may have been made to the software and hardware it describes since then. ADInstruments NZ Limited reserves the right to alter specifications as required. Late-breaking information may be supplied separately.

Trademarks of ADInstruments

PowerLab®, LabChart® and ADInstruments® are registered trademarks of ADInstruments NZ Limited.

Other Trademarks

Apple, Mac and Macintosh are registered trademarks of Apple Computer, Inc.

Windows, Windows 7, Windows 8, Windows 10 and Windows Vista are either registered trademarks or trademarks of Microsoft Corporation.

All other trademarks are the property of their respective owners.

Document Number: U-BA/OG-01 Rev D, Date of issue: 04/22

Copyright © ADInstruments NZ Limited, 2022. All rights reserved. PowerLab, LabChart and ADInstruments are registered trademarks of ADInstruments NZ Limited. Windows 8, Windows 7, Windows 10, Windows Vista and .NET Framework are trademarks of Microsoft Corporation. Apple, the Apple logo, MacOS, and Macintosh are trademarks of Apple Computer Inc. registered in the U.S. and other countries. Acrobat and Adobe are registered trademarks of Adobe Systems Incorporated. Igor is a trademark of Wavemetrics Inc. MATLAB is a registered trademark of The MathWorks Inc. Grass is a trademark of Astro-Med Inc. All other trademarks are the property of their respective owners.

Web: www.adinstruments.com

Manufactured in Australia by:
ADInstruments (Sydney) Pty. Ltd.,
13/22 Lexington Drive
Bella Vista 2153 New South Wales

Technical Support: support.au@adinstruments.com



Chapter 1

Safety Notes

Statement of Intended Use

All products manufactured by ADInstruments are intended for use in teaching and research applications and environments only. ADInstruments products are NOT intended to be used as medical devices or in medical environments. That is, no product supplied by ADInstruments is intended to be used to diagnose, treat or monitor a subject. Furthermore no product is intended for the prevention, curing or alleviation of disease, injury or handicap. ADInstruments products are intended to be installed, used and operated under the supervision of an appropriately qualified life-science researcher. The typical usage environment is a research or teaching lab or hospital. ADInstruments equipment is not intended for use in domestic environments.

Where a product meets IEC 60601-1 it is under the principle that:

- this is a more rigorous standard than other standards that could be chosen.
- it provides a high safety level for subjects and operators.

The choice to meet IEC 60601-1 is in no way to be interpreted to mean that a product:

- is a medical device,
- may be interpreted as a medical device, or
- is safe to be used as a medical device.

Safety and Quality Standards

When used with ADInstruments isolated front-ends, PowerLab systems are safe for connection to subjects. The FE231 Bio Amp, FE232 Dual Bio Amp and FE234/FE238 Quad/Octal Bio Amps front-ends conform to international safety requirements. Specifically these are IEC60601-1 and its addenda (Safety Standards, page 3) and various harmonized standards worldwide (CSA601.1 in Canada and AS/NZS 3200.1 in Australia and New Zealand).

In accordance with European standards they also comply with the electromagnetic compatibility requirements under IEC60601-1-2, which ensures compliance with the EMC directive.

Quality Management System ISO 9001:2008

ADInstruments manufactures products under a quality system certified as complying with ISO 9001:2008 by an accredited certification body.

Regulatory Symbols

Amplifiers and signal-conditioners manufactured by ADInstruments that are designed for direct connection to humans and animals are tested to IEC60601-1:2012 (including amendments 1 and 2), and carry one or more of the safety symbols below. These symbols appear next to those inputs and output connectors that can be directly connected to human subjects.



BF (body protected) symbol. This means that the input connectors are suitable for connection to humans and animals provided there is no direct electrical connection to the heart.



Warning symbol. The exclamation mark inside a triangle means that the supplied documentation must be consulted for operating, cautionary or safety information before using the device.



CE Mark. All front-end amplifiers and PowerLab systems carry the CE mark and meet the appropriate EU directives.



UKCA mark. All front-end amplifiers and PowerLab systems carry the UKCA mark and meet the appropriate UK directives.



Refer to booklet symbol. This symbol specifies that the user needs to refer to the Instruction manual or the booklet associated with the device.



Date of Manufacture/ Manufacturer's name symbol. This symbol indicates the date of manufacture of the device and the name of the manufacturer



WEEE directive symbol. Unwanted equipment bearing the Waste Electrical and Electronic Equipment (WEEE) Directive symbol requires separate waste collection. (See disposal section at the end of this chapter)

Further information is available on request.

Safety Standards

IEC Standard - International Standard - Medical Electrical Equipment

IEC 60601-1-1:2000 Safety requirements for medical electrical systems

IEC 60601-1:2012 + A1 General requirements for safety

General Safety Instructions

To achieve the optimal degree of subject and operator safety, consideration should be given to the following guidelines when setting up a PowerLab system either as stand-alone equipment or when using PowerLab equipment in conjunction with other equipment. Failure to do so may compromise the inherent safety measures designed into PowerLab equipment. ADInstruments front-ends are only suitable for operation with ADInstruments PowerLabs. Front-ends are suitable for use with any S/, SP/, /20, /25, /30 and /35 series and 15T PowerLabs (FE234 and FE238 only suitable for use with 35 series PowerLabs). Note that compliance with IEC60601-1 can only be achieved when front-ends are used with a /35 series Powerlab.

The following guidelines are based on principles outlined in the international safety standard IEC 60601-1: *General requirements for safety – Collateral standard: Safety requirements for medical systems*. Reference to this standard is required when setting up a system for human connection. The user is responsible for ensuring any particular configuration of equipment complies with IEC60601-1-1. Guidance on compliance with this standard is provided in the following sections.

PowerLab systems (and many other devices) require the connection of a personal computer for operation. This personal computer should be certified as complying with IEC 60950 and should be located outside a 1.8 m radius from the subject (so that the subject cannot touch it while connected to the system). Within this 1.8 m radius, only equipment complying with IEC 60601-1 should be present. Connecting a system in this way obviates the provision of additional safety measures and the measurement of leakage currents.

Accompanying documents for each piece of equipment in the system should be thoroughly examined prior to connection of the system.

While it is not possible to cover all arrangements of equipment in a system, some general guidelines for safe use of the equipment are presented below:

- Any electrical equipment which is located within the SUBJECT AREA should be approved to IEC 60601-1.
- Only connect those parts of equipment that are marked as an APPLIED PART to the subject. APPLIED PARTS may be recognized by the BF symbol which appears in the Safety Symbols section of these Safety Notes.
- Never connect parts which are marked as an APPLIED PART to those which are not marked as APPLIED PARTS.

-
- Do not touch the subject to which the PowerLab (or its peripherals) is connected at the same time as making contact with parts of the PowerLab (or its peripherals) that are not intended for contact to the subject.
 - Cleaning and sterilization of equipment should be performed in accordance with manufacturer's instructions. The isolation barrier may be compromised if manufacturer's cleaning instructions are not followed.
 - The ambient environment (such as the temperature and relative humidity) of the system should be kept within the manufacturer's specified range or the isolation barrier may be compromised.
 - The entry of liquids into equipment may also compromise the isolation barrier. If spillage occurs, the manufacturer of the affected equipment should be contacted before using the equipment.
 - Many electrical systems (particularly those in metal enclosures) depend upon the presence of a protective earth for electrical safety. This is generally provided from the power outlet through a power cord, but may also be supplied as a dedicated safety earth conductor. Power cords should never be modified so as to remove the earth connection. The integrity of the protective earth connection between each piece of equipment and the protective earth should be verified regularly by qualified personnel.
 - Avoid using multiple portable socket-outlets (such as power boards) where possible as they provide an inherently less safe environment with respect to electrical hazards. Individual connection of each piece of equipment to fixed mains socket-outlets is the preferred means of connection.

If multiple portable socket outlets are used, they are subject to the following constraints:

- They shall not be placed on the floor.
- Additional multiple portable socket outlets or extension cords shall not be connected to the system.
- They shall only be used for supplying power to equipment which is intended to form part of the system.

Bio Amp Safety Instructions

The Bio Amp inputs displaying any of the safety symbols are electrically isolated from the mains supply in order to prevent current flow that may otherwise result in injury to the subject. Several points must be observed for safe operation of the Bio Amp:

- All Bio Amp front-ends (except for the FE234 Quad and FE238 Octal Bio Amps) and all PowerLab units with a built-in Bio Amp are supplied with a 3-lead or 5-lead Bio Amp subject cable and lead wire system. The FE234 Quad and FE238 Octal Bio Amps are supplied with unshielded lead wires (1.8 m). Bio Amps are only safe for human connection if used with the supplied subject cable and lead wires.
- Bio Amp front-ends are NOT defibrillator-protected. Using the Bio Amp to record signals during defibrillator discharges may damage the input stages of the amplifiers. This may result in a safety hazard.
- Never use damaged Bio Amp cables or leads. Damaged cables and leads must always be replaced before any connection to humans is made.

Earthing and Ground Loop Noise

The prime function of earthing is safety, that is, protection against fatal electrocution. Safety concerns should always override concerns about signal quality. Secondary functions of earthing are to provide a reference potential for the electrical equipment and to mitigate against interference.

The earthing (grounding) stud provided on the back panel of the PowerLab is a potential equalization post and is compatible with the DIN 42801 standard. It is directly connected to the earth pin of the power socket and the PowerLab chassis. The earthing stud can be used where other electronic equipment is connected to the PowerLab, and where conductive shields are used to reduce radiative electrical pick-up. Connection to the stud provides a common earth for all linked devices and shields, to reduce ground-loops.

The earthing stud can also be used where a suitable ground connection is not provided with the mains supply by connecting the stud to an earthed metal infrastructure, such as a metal stake driven into the ground, or metal water piping. This may also be required in laboratories where safety standards require additional grounding protection when equipment is connected to human subjects. Always observe the relevant safety standards and instructions.

Note that electromagnetically-induced interference in the recorded signal can be reduced by minimizing the loop area of signal cables, for example by twisting them together, or by moving power supplies away from sensitive equipment to reduce the inductive pick-up of mains frequency fields. Please consult a good text for further discussion of noise reduction.

Cleaning and Sterilization

ADInstruments products may be wiped down with a lint free cloth moistened with industrial methylated spirit. Refer to the manufacturer's guidelines or the Data Card supplied with transducers and accessories for specific cleaning and sterilizing instructions.

Inspection and Maintenance

PowerLab systems and ADInstruments front-ends are all maintenance-free and do not require periodic calibration or adjustment to ensure safe operation. Internal diagnostic software performs system checks during power up and will report errors if a significant problem is found. There is no need to open the instrument for inspection or maintenance, and doing so within the warranty period will void the warranty.

Your PowerLab system can be periodically checked for basic safety by using an appropriate safety testing device. Tests such as earth leakage, earth bond, insulation resistance, subject leakage and auxiliary currents and power cable integrity can all be performed on the PowerLab system without having to remove the covers. Follow the instructions for the testing device if performing such tests. If the PowerLab system is found not to comply with such testing you should contact your PowerLab representative to arrange for the equipment to be checked and serviced.



WEEE Directive
symbol

Environment

Electronic components are susceptible to corrosive substances and atmospheres, and must be kept away from laboratory chemicals.

Disposal

- Forward to recycling center or return to manufacturer.
- Unwanted equipment bearing the Waste Electrical and Electronic Equipment (WEEE) Directive symbol requires separate waste collection. For a product labeled with this symbol, either forward to a recycling center or contact your nearest ADInstruments representative for methods of disposal at the end of its working life.



Chapter 2

Overview

The PowerLab system consists of a recording unit and application programs that run on the computer to which the unit is connected. It provides an integrated system of hardware and software designed to record, display, and analyze experimental data.



Front-ends are ancillary devices that connect to the PowerLab recording unit to extend the system's capabilities. They provide additional signal conditioning, and other features, and extend the types of experiments that you can conduct and the data you can record.

All ADInstruments front-ends are designed to be operated under full software control. No knobs, dials, or switches are needed, although some may be provided for reasons of convenience or safety.

Introduction

The PowerLab controls front-ends through an expansion connector called the I²C (eye-squared-sea) bus. This makes it very easy to add front-ends to the system or to transfer them between PowerLabs. Many front-ends can be added to the system by connecting the I²C sockets in a simple daisy-chain structure. The PowerLab provides control and low-voltage power to front-ends through the I²C bus so, in general, no separate power supply is required.

In addition, each front-end requires a separate connection to one or more analog input channel(s) of the PowerLab. External signals are acquired through the PowerLab analog inputs and amplified before being digitized by the PowerLab. The digitized signal is transmitted to the computer using a fast USB connection. ADInstruments software applications LabChart, LabTutor, LabStation and Lt receive, display, and record the data and your analysis to the computer's hard disk.

Front-ends are automatically recognized by the PowerLab system. Once connected, the features of the front-end are combined with the appropriate features of the PowerLab (for example, range and filtering options) and are presented as a single set of software controls.

Note: The Stimulator front-ends differ from other front-ends in two respects:

1. Since they need to produce a reasonably high voltage and current, the Stimulator front-ends require a power supply in addition to the power provided by the I²C bus.
2. As they produce voltage output for stimulation, they are connected to a positive analog output socket of the PowerLab as a source for timing and producing pulses.

A variety of accessory products are available with ADInstruments Front-ends, such as transducers, signal cables and recording electrodes. Some of these are listed in the Getting Started with Front-end Signal Conditioners booklet, supplied with your Front-end. For more details see: <http://www.adinstruments.com/> or contact your local ADInstruments representative.

Checking the Front-end

Before connecting the front-end to anything, check it carefully for signs of physical damage.

1. Check that there are no obvious signs of damage to the outside of the front-end casing.
2. Check that there is no obvious sign of internal damage, such as rattling. Pick up the front-end, tilt it gently from side to side, and listen for anything that appears to be loose.

If you have found a problem, contact your authorized ADInstruments representative immediately and describe the problem. Arrangements can be made to replace or repair the front-end.

Connecting to the PowerLab


To connect a front-end to the PowerLab, first ensure that the PowerLab is turned off. Failure to do this may damage the PowerLab, the front-end, or both.

The BNC cable from the front-end signal output must connect to an analog input on the PowerLab. If you have an older PowerLab that has differential (rather than single-ended) inputs, the front-end must connect to a *positive* input.

Single Front-ends

Connect the I²C output of the PowerLab to the I²C input of the front-end using the I²C cable provided. Figure 2–1 shows how to connect up a single front-end to your recording unit.

Figure 2–1
Connecting a front-end to the PowerLab: a PowerLab has only one I²C output, and each front-end has one I²C output and one I²C input



Check that the connectors for the I²C bus are screwed in firmly. Check the BNC cable for firm connections as well. Loose connectors can cause erratic front-end behavior, or may cause the front-end to fail to work at all.

The Signal Output Socket


The BNC socket labelled Signal Output on the back panel of the front-end provides the signal output to connect to an analog input socket on the front of the PowerLab. A BNC-to-BNC cable is supplied for this connection. If necessary, use a BNC to DIN smart adapter [MLAC22] to connect the BNC cable to your PowerLab's input.

Note: If you have an older PowerLab with differential (rather than single-ended) inputs, the BNC cable must connect to a *positive* analog input on the PowerLab.

Multiple Front-ends

Multiple separate front-ends can be connected up to a PowerLab. The initial front-end should be connected with the I²C cable as in Figure 2-1. The remainder are daisy-chained via I²C cables, connecting the I²C output of the last connected front-end to the I²C input of the front-end to be added (Figure 2-2).

Figure 2-2
Connecting multiple front-ends to the PowerLab (two single front-ends shown for simplicity)



The number of normal front-ends that can be connected to a PowerLab depends on the number of analog input channels on the PowerLab. Each BNC cable from a front-end should be connected to one analog input channel on the PowerLab, for example, Input 1 on a /30 or /35 series PowerLab.

Note: Only one Stimulator front-end such as a Stimulus Isolator can be connected to the positive output of the PowerLab.

Special Cases

Some front-ends have their own specific connection requirements. Please refer to the individual chapter for each front-end in this guide.

Connecting Stimulator Front-Ends

The PowerLab analog outputs provide a variable, computer-controlled voltage output that can be used with LabChart, LabTutor, LabStation or Lt to connect a Stimulator front-end, or to stimulate directly, or to control a peripheral device. A voltage output is generated by the PowerLab and delivered via the BNC output sockets, giving positive, negative, differential, or independent stimuli, depending on the PowerLab used and the software settings.

The /20, /25, and /26 series PowerLabs have analog outputs labeled + and -. In contrast, the SP, ST, /30 and /35 series PowerLabs have the outputs labeled Output 1 and Output 2.

For the /20, /25 and /26 series PowerLabs:

The negative (–) output is the complement of the positive (+) output, so the stimuli from the two outputs are mirror images. If one output gives a positive voltage, the other gives a negative one, and the two together give a differential voltage. One Stimulator front-end such as a Stimulus Isolator or Stimulator HC can be connected to the positive output of these PowerLabs.

Note: If you connect the Stimulator HC to a PowerLab that has an in-built Isolated Stimulator, such as a PowerLab 26T, only the external, connected stimulator is used.

For /SP, /ST, /30 and /35 series PowerLabs:

Output 1 and Output 2 can function independently. However, only one Stimulator front-end such as a Stimulus Isolator or Stimulator HC can be connected to the positive output (Output 1) of these PowerLabs. With a Stimulator front-end connected, the second output (Output 2) can function independently, and a second tab appears in the Stimulator dialog in LabChart 7 for Windows. Therefore Output 2 remains available for other uses, such as creating analog waveforms and triggering other systems.

Maximum Number of Front-Ends

The I²C bus can control a maximum of sixteen front-ends. Therefore, if you are using a PowerLab 16/30, which has sixteen input channels, you can record from sixteen single channel front-ends.

Using ADInstruments Programs

Front-ends are designed for use with PowerLabs and ADInstruments programs such as LabChart, LabTutor, LabStation and Lt. The functions of the front-end are combined with those of the PowerLab, and are presented as a single set of software controls in the ADInstruments program. Depending on the front-end(s) connected, front-end-specific dialogs replace the Input Amplifier dialogs or the Stimulator dialog.

The **LabChart Help** detail the Input Amplifier and Stimulator dialogs, and explain relevant terms and concepts, but they do not cover front-end-specific features. These features are described in detail in the following chapters for each front-end.

Front-end Drivers

A device driver is a piece of software that allows the computer's operating system and other software to interact with a hardware device. ADInstruments applications like LabChart communicate with a front-end via an appropriate front-end driver. These drivers are automatically set up on the computer when ADInstruments applications are installed, and their operation is usually invisible to the user.

However, under certain circumstances you may receive an error message during the startup of LabChart indicating that there is a problem with the front-end driver. Subsequently, the front-end will not function. This is invariably caused by the absence or incompatibility of a driver required for communication with the front-end due to an old version of the software being run. The problem can be remedied simply by reinstalling

and rerunning a current version of the software, which will include the latest front-end drivers.

The Front-end Self-test

Once the front-end is properly connected to the PowerLab, and the proper software is installed on the computer, a quick check can be performed on the front-end. To perform the self-test:

- Turn on the PowerLab and check that it is working properly, as described in the owner's guide that was supplied with it.
- Once the PowerLab is ready, start LabChart, LabTutor, LabStation or Lt.
- While the program is starting, watch the Status indicator on the front-end's front panel. During initialization, you should see the indicator flash briefly and then remain lit.

If the indicator lights correctly, the front-end has been found by the PowerLab and is working properly. If the indicator doesn't light, check your cable connections and repeat the start-up procedure.


Software Behavior

When a front-end is connected to a PowerLab and the ADInstruments software is successfully installed, the **Input Amplifier...** menu command from the Channel Function pop-up menu in LabChart should be replaced by the **<Front-end>...** menu command.

For example, with a Bio Amp front-end connected, **Bio Amp...** should appear in the Channel function pop-up menu.

Figure 2-3

Channel Function pop-up menu in LabChart with the Bio Amp front-end connected



If the application fails to find a front-end attached to a channel, the normal **Input Amplifier...** command or button remains. If you were expecting a connected front-end, you should close the program, turn everything off, check the connections, restart the PowerLab and then relaunch LabChart, LabTutor or the Kuraloud Desktop App.

Preventing Problems

Several problems can arise when using the PowerLab system for recording biological signals. It is important to understand the types of problems that can occur, how they manifest themselves, and what can be done to remove them or to minimize their effect. These are usually problems of technique, and should be addressed before you set up your equipment.

Aliasing

Recordings of periodic waveforms that have been undersampled may have misleading shapes and may also have artifacts introduced by aliasing. Aliasing occurs when a regular signal is digitized at too low a sampling rate, causing the false appearance of lower frequency signals. An analogy to aliasing can be seen in old films: spoked wagon wheels may appear to stop, rotate too slowly or even go backwards when their rate of rotation matches the film frame speed – this is obviously not an accurate record.

The Nyquist-Shannon sampling theorem states that the minimum sampling rate (f_s) to accurately describe an analog signal must be at least twice the highest frequency in the original signal. Therefore, the signal must not contain components greater or equal to $f_s/2$. The term $f_s/2$ is known as the Nyquist frequency (f_n) or the ‘folding frequency’ because frequencies greater than or equal to f_n fold down to lower frequencies about the axis of f_n .

When aliasing of noise or signals is seen, or even suspected, the first action you should take is to increase the sampling rate. The highest available sampling rates are 100k /s or 200k /s, depending on your PowerLab. To view the frequencies present in your recorded signal open the Spectrum window in LabChart. For more information about Spectrum, see the LabChart Help Center.

If unwanted high-frequency components are present in the sampled signal, you will achieve better results by using a low-pass filter to remove them. The best kind of filter for this purpose is the Anti-alias filter option available in the front-end-specific **Input Amplifier...** dialog. This is a special low-pass filter that is configured to automatically remove all signals that could alias; i.e., those whose frequency is greater or equal to half the sampling rate.

For certain PowerLabs, the Anti-alias filter option is not available. Therefore you should select an appropriate low-pass filter to remove any unwanted signals (or noise) occurring at frequencies greater or equal to half the sampling rate.

Frequency Distortion

Frequency distortion will occur if the bandwidth of your recording is made smaller than the bandwidth of the incoming signal. For example, if an ECG was measured with a sampling rate of 100 samples per second (100 Hz) and the Bio Amp had a low-pass filter applied at 50 Hz, the fast-changing sections of the waveform (the QRS complex) may appear smaller and ‘blunted’, while the slower T-wave sections remain relatively unchanged. This overall effect is called frequency distortion.

It can be eliminated by increasing the frequency cut-off of the low-pass filter in the front-end-specific **Input Amplifier...** dialog to obtain an undistorted waveform.

Similarly, if the high-pass filter was set too high, the amplitude of the T-wave sections may be reduced. The **Input Amplifier...** dialog allows you to examine ECGs and similar slowly changing waveforms to fine-tune filter settings before recording.

Saturation

Saturation occurs when the range is set too low for the signal being measured (the amplification, or gain, is too high). As the signal amplitude exceeds the allocated range, the recorded waveform appears as if part of the waveform had been cut off, an effect referred to as clipping.

Clipping can also be caused by excessive baseline offset: the offset effectively moves the whole waveform positively or negatively to an extent that causes all or part of it to be clipped. This problem is overcome by selecting a higher range from the Range menu in the front-end-specific **Input Amplifier...** dialog. In the case of excessive baseline offset, you may wish to apply a high-pass filter with a higher frequency cut-off.

Ground Loops

Ground loops occur when multiple connected pieces of recording equipment are connected to mains power grounds. For safety reasons, *all* electrical equipment should have a proper connection to the mains power grounds, or to a primary earth connection in situations where a mains ground connection is not available. Connecting linked electrical equipment to a common earth connection (equipotential connection point) – such as the earthing (grounding) stud provided on the rear of all PowerLabs – can prevent ground loops.

The electric fields generated by power lines can introduce interference at the line frequency into the recorded signal. Electromagnetic fields from other sources can also cause interference: fluorescent tubes, apparatus with large transformers, computers, laptop batteries, network cables, x-ray machines, microwave ovens, electron microscopes, even cyclic air conditioning.

Reasonable care in the arrangement of equipment to minimize the ground loop area, together with proper shielding, can reduce electrical frequency interference. For example, use shielded cables, keep recording leads as short as possible, and try twisting recording leads together. For sensitive measurements, it may be necessary to place the subject (the biological source) in a Faraday cage.

Interference should first be minimized, and then you can turn on the Mains filter in the front-end-specific **Input Amplifier...** dialog.

Mains filter

The Mains filter (/20, /25, /30, /35 and 26T PowerLabs) allows you to filter out interference at the mains frequency (typically 50 or 60 Hz). The mains filter is an adaptive filter which tracks the input signal over approximately 1 second. A template of mains-frequency signal present in the input is computed from the signal. The width of the template is the mains power period (typically 16.6 or 20 ms) as determined from zero-crossings of

the mains power. The filtered signal is obtained by subtracting the template from the incoming signal.

In comparison with a conventional notch filter, this method produces little waveform distortion. It attenuates harmonics of the mains frequency as well as the 50 or 60 Hz fundamental and therefore effectively removes non-sinusoidal interference, such as that commonly caused by fluorescent lights.

The filter should not be used when:

- the interference changes rapidly. The filter takes about 1 second to adapt to the present level. If interference is present and then is suddenly removed, interference in the filtered signal will temporarily worsen.
- your signal contains exact factors or harmonics of frequencies close to the mains frequencies, for example, a 30 Hz signal with 60 Hz mains frequency.
- your signal is already free from interference. If the signal-to-noise ratio is greater than about 64 the mains filter introduces more noise than it removes.
- you are recording at close to maximum sampling rates. The mains filter uses some of the PowerLab's processing power and therefore reduces the maximum rate at which you can sample.

Electrode Contact


Occasionally one of the lead wires connecting the subject to the front-end may become disconnected, or an electrode contact may become poor. If this should happen, relatively high voltages (potentials) can be induced in the open wire by electric fields generated by power lines or other sources close to the front-end or the subject. Such induced potentials will result in a constant amplitude disturbance in the recorded waveform at the power line frequency (50 or 60 Hz), and loss of the desired signal. If the problem is a recurring one, one of the leads may be faulty. Check connections and replace faulty leads, if necessary.

Motion Artifacts

A common source of artifacts when recording biological signals is due to motion of the subject or equipment. Often applying a high-pass filter can help to remove slowly changing components in a recorded signal.

- Muscular activity generates its own electrical signals, which may be recorded along with an ECG, say, depending on the location of the electrodes.
- If an electrode is not firmly attached, impedance (and hence the recorded signal) may vary as the contact area changes shape owing to movement.
- Movement of patient cables, particularly bending or rubbing together (triboelectric effects) may generate artifacts in a signal.
- Subject respiration can also generate a signal; breathing can result in a slowly changing baseline corresponding to inspiration and expiration.

If the subject is liable to move during recording, then special care needs to be taken when attaching the electrodes and securing the patient leads. Make sure the skin is cleaned and lightly abraded before attaching the electrodes.





Chapter 3

Bio Amp

This chapter provides an overview of the Bio Amp [FE231], Dual Bio Amp [FE232], Quad Bio Amp [FE234] and Octal Bio Amp [FE238]. The Bio Amp allows the PowerLab system to record biological signals, such as ECGs, EMGs, and EEGs from humans or animals.



Bio Amp Safety Instructions

WARNING:
Refer to Intended
Use statement
on page 1 of this
Owner's Guide
before use.

The Bio Amp inputs displaying safety symbols are electrically isolated from the mains supply in order to prevent current flow that may otherwise result in injury to the subject. Several points must be observed for safe operation of the Bio Amp:

- The Bio Amp [FE231], Dual Bio Amp [FE232] and PowerLab units with a built-in Bio Amp are supplied with a 3-lead or 5-lead Bio Amp subject cable and lead wire system. The Octal Bio Amp [FE238] and Quad Bio Amp [FE234] are supplied with unshielded lead wires (1.8m). Bio Amps are only safe for human connection if used with the supplied subject cable and lead wires.
- All Bio Amp front-ends and PowerLab units with a built-in Bio Amp are not defibrillator-protected. Using the Bio Amp to record signals during defibrillator discharges may damage the input stages of the amplifiers. This may result in a safety hazard.

- Never use damaged Bio Amp cables or leads. Damaged cables and leads must always be replaced before any connection to humans is made.


The Bio Amp front-ends

All Bio Amps consist of electrically isolated, differential input AC/DC amplifiers with common isolated ground connections. The number of amplifiers corresponds to the number of Bio Amp channels. For example, the Dual Bio Amp has two amplifiers with a shared ground connector and input socket, and the Octal Bio Amp has eight differential amplifiers with a shared ground connector and separate inputs. For multiple biological recordings (more than two) from a single subject, it is recommended that you use a multi-channel Bio Amp rather than “daisy-chaining” multiple Bio Amps.

The Front Panel

The front panel of the Single Bio Amp has three input connectors and one indicator light.

Figure 3-1
The front views of
the Bio Amp, Dual
Bio Amp, Quad
Bio Amp and
Octal Bio Amp



WARNING:
When used
in ambient
temperatures
of 38 degrees C
and above, do
not touch the Bio
Amp enclosure
for more than
a minute
continuously.

The front panel of a Dual Bio Amp has a single input, two indicator lights and an audio output connector.

The front panel of the Quad Bio Amp has nine single-ended 1.5mm sockets consisting of four red connectors, four black connectors and a single green connector.

The front panel of the Octal Bio Amp has seventeen single-ended 1.5 mm sockets consisting of eight red connectors, eight black connectors and a single green connector.

The Input Socket

Connections are made to the Bio Amp are made using the three shrouded 1.5 mm male pin sockets on the front panel. A separate socket is provided for each of the positive (+), negative (-) and Ground/Reference (Ref) cables.

Connections are made to the Dual Bio Amp using the six-pin socket on the front panel. The socket is physically and electrically isolated from the low-voltage mains-supply circuitry of the PowerLab, and the input connections are isolated internally, by isolation circuitry. The socket is of a sort commonly used with ECG-type cables and leads, such as the Bio Amp cable and leads with which your Bio Amp is supplied (Tronomed D-1340 or Tronomed D-1540).

Connections are made to the Quad/Octal Bio Amp using individual lead wires (supplied with every unit) and a common ground, rather than using a single multi-pin socket. Each of the Bio Amp inputs has one red and one black connector and an indicator light. The green connector is a shared ground connection across all eight inputs in Octal and four in Quad.

The Status Indicator

The status indicator light of a Single or Dual Bio Amp is located at the bottom left of the front panel. The status indicator lights of an Octal and Quad Bio Amp are located above the pair of red and black connectors of each input on the front panel.

When an ADInstruments application such as LabChart starts, the status indicators should flash briefly and then remain green, indicating that the program has found the front-end, checked and selected it, and is ready to use it. If a status indicator does not turn on and stay on when the application starts, it is most likely that the front-end is not connected properly.

The Back Panel

The back panel of the Bio Amp provides all the sockets required to connect the front-end to the PowerLab and to other front-ends.

Audio Out Socket

Bio Amps have audio monitor outputs on the rear panel that can be used with standard headphones or externally powered speakers. The 3.5mm headphone socket provides sound output from a software-selected data channel. The audio output may be of use when monitoring nerve firings to control the placement of electrodes, for instance.


I²C Input and Output

The Bio-Amp communicates with the PowerLab via the 9 Pin I²C Input connector. This connector can either be connected directly to a PowerLab or to the 'I²C Output' connector of another Front-End. Up to 16 Front Ends can be 'daisy chained' in this fashion. More details can be found in the 'Connecting to the PowerLab' section in Chapter 2. The I²C socket allows multiple front-ends to be used independently with one PowerLab. Power and control signals to connected front-ends come from the PowerLab.

Figure 3-2

The rear views of the Bio Amp, Dual Bio Amp, Quad Bio Amp and Octal Bio Amp

WARNING:
All Bio Amplifiers have only been assessed as compliant with IEC60601-1 safety standard when used with a 35 series PowerLab



Analog Output Sockets

The BNC sockets on the back panel of the Bio Amps provide the signal outputs to connected analog input sockets of the PowerLab. The sockets are labeled Signal Output on a Single Bio Amp and Output 1 to Output n on a multi-channel Bio Amp. A BNC-to-BNC cable is supplied for each connection. It is recommended that, to avoid confusion, you match the Bio Amp output number with the corresponding PowerLab input connector.

Audio Out Socket

All Bio Amps have an audio monitor output on the back panel that can be used with a wide range of headphones or externally powered speakers. The 3.5 mm stereo socket provides mono sound. The audio output may be of use when monitoring nerve firings to control the placement of electrodes, for instance.

WARNING:

As with all audio devices, there is a risk of temporary hearing damage if used carelessly. It is recommended headphones are attached to the Neuro Amp or Bio Amp audio socket prior to attaching the headphones to ears.



Applied parts

Single/Dual Bio Amp:

- MLA2540 (Dual) Bio Amp Cable
- MLA2505 Lead wires
- MLA1212 Micro-Hook Electrodes (Single)

Octal/Quad Bio Amp:

- MLA0310 Lead wires, unshielded

Other Supplied Accessories

Single/Dual Bio Amp:

- MLAC01 BNC-BNC cables
- MLAC02 DB9M- DB9M cable
- MLA2350 (Single) Bio Amp Cable

Octal/Quad Bio Amp:


- MLAC27 Quad BNC-BNC cables
- MLAC02 DB9M- DB9M cable

The Bio Amp Cable

Connections should be made to the Bio Amp and Dual Bio Amp inputs using the supplied Bio Amp cables and leads. The Dual Bio Amp cable plugs into the six-pin input socket on the Dual Bio front panel: a notch in the plug ensures that polarity is correct. The lead wires on the Bio Amp cable plug into the input sockets on the Bio Amp front panel. Only the supplied Bio Amp cable and leads should be used. Other cables may not meet safety requirements.

Figure 3-3

Bio Amp cable yokes, with leads attached: 3 leads for the Bio Amp, and 5 leads for the Dual Bio Amp



The single Bio Amp is supplied with a 3-lead Bio Amp cable with lead wires. The Dual Bio Amp is supplied with a 5-lead Bio Amp cable and separate lead wires; it uses a shared ground signal for its Bio Amp inputs. The supplied cables are of the type used for ECG or EMG studies: Tronomed D-1340 or Tronomed D-1540 cables respectively. The cable has a yoke with three or five sockets for the lead wires. Note that the active pins of a shielded cable are those closest to the label.

The lead wires supplied are of the type used for ECG studies. They click into place in the cable yoke, and have snap connectors at the other end to connect to typical ECG electrodes. The leads are color-coded for identification. The labels on the Bio Amp cable also have color spots to help sort out which cables connect where and what they are measuring. (The colors are arbitrary, since the PowerLab system is for general-purpose recording.)


Connections to the Octal and Quad Bio Amp inputs can be made with the supplied lead wires. The input sockets are for unshielded, single-ended lead wires with a 2 mm inside diameter with 1.5 mm pins.

ADInstruments supplies other types of lead wires that connect to the Bio Amp cable yoke. Also available are disposable and reusable electrodes, electrode cream (for reusable electrodes), and abrasive gel, for lightly abrading the skin before the electrodes are attached.

The Bio Amp Cable Input

Connections are made to the Bio Amp using the three shrouded 1.5 mm male pin sockets on the front panel. The sockets are of a sort commonly used with life science connection leads, and their arrangement is shown below.

Two types of leads are provided: one set of 3 individual leads terminating in alligator clips (commonly used for animal or in vitro recordings) and one shielded Bio Amp cable with detachable snap-on leads (commonly used for human recordings).



The sockets provide two pins for a differential input signal (+ or red, - or black), and an input for an isolated Ground/reference (Ref or green). The Ground/reference electrode should always be attached to the subject under investigation (or to the recording preparation), providing a zero reference for the differential amplifier. Then the active (+) and reference (-) electrodes should be securely attached to the subject or recording preparation, for example, ECG recording leads on opposite sides of the chest, or needle electrodes at different sites over a nerve.

The Bio Amp cable for the Dual Bio Amp, plugs into the six-pin input socket on the front panel of the Dual Bio Amp. A notch in the plug ensures that polarity is correct. Only the


Figure 3-4
The pin assignments for the single Bio Amp input connector

supplied Bio Amp cable and lead wires should be used as other cables may not meet safety requirements.

The Dual Bio Amp has one common connector for two Bio Amp channels, nominally channels 1 and 2. The biological amplifiers both have differential inputs, a shared Ground/reference, and connection to the cable's shield. The entire connector is physically and electrically isolated to ensure subject safety.

Figure 3-5


The pin assignments for the Dual Bio Amp input connector



The Octal Bio Amp has a pair of connectors for each of the eight Bio Amp inputs, and a single connector for the shared Ground/reference. The Quad has a pair of connectors for each the four Bio Amp inputs, and a single connector for the shared Ground/reference. The connectors are physically and electrically isolated to ensure subject safety. The Octal and Quad Bio Amps should be used with the supplied 1.8m unshielded lead wires.

Figure 3-6

The connectors for each input of the Octal/Quad Bio Amp



The Common lead must always be connected to the subject to ensure that the common-mode signal seen by the inputs is not unduly large. An unduly large common mode signal may lead to poor signal quality.

Shutdown Procedure: After use, the Bio Amp should be shut down by quitting LabChart and powering off the PowerLab

Using LabChart

Once the Bio Amp is connected, turn the PowerLab on and launch LabChart. When a Bio Amp is properly connected to the PowerLab, the **Input Amplifier...** menu command is replaced by **Bio Amp...** for the input channel to which it is connected.

If LabChart fails to find a front-end connected, the normal text remains. If you were expecting a connected front-end and see the normal text instead, you should quit

the application, turn the PowerLab off and check the connections. Then restart the PowerLab and re-launch LabChart.


The documentation for LabChart does not cover front-end-specific features. These features are described in detail here for LabChart.

The Bio Amp dialog

The Bio Amp dialog (Figure 4–8 and Figure 4–9) allows software control of the combined input amplifiers and filters in the PowerLab and Bio Amp. The signal present at a channel's input is displayed in the preview area so that you can see the effects of changes in settings. Change settings in the dialog, then click the OK button to apply them.

The Bio Amp dialog appears when you choose **Bio Amp...** from a Channel Drop Down menu. To set up many channels quickly, open the **Setup > Channel Settings...** dialog. Here you can view all the channels that are turned on, and you can turn off any unnecessary channels. Clicking on **Bio Amp...** in the Input Settings column of the Channel Settings dialog will also open the Bio Amp dialog.

Figure 3–7
The Bio Amp
dialog in LabChart
for Windows for
the Dual
and Octal/Quad
Bio Amps



Signal Display

The input signal is displayed so you can see the effect of changing the settings — no data are recorded while setting things up. The instantaneous signal value is shown above the display area.

Setting the Range

The Range pop-up menu lets you select the input range of the channel (combined range of the PowerLab and Bio Amp). Note that range is inversely related to gain (e.g, 100 mV range is the lowest gain and sensitivity setting). Changing the range in the Bio Amp dialog is equivalent to changing it using the arrow in the top right hand corner of the channel in the Chart View Window. Changing the hardware range setting should not be confusing with adjusting the y axis scaling using the + and - buttons on the left of each channel. The default setting (if you have not loaded a settings file) is 100 mV and the ranges go down to 100 μ V in 10 steps.

Filtering

The Bio Amp has low-pass, high-pass, and mains filter circuitry that can be adjusted to suit the application. The mains filter removes excessive mains-frequency interference. The high-pass and low-pass filters provide bandwidth limiting of low-frequency and high-frequency signals respectively. Note that the settings for one filter type may restrict the possible settings for the other.

High-Pass Filtering. The High Pass pop-up menu gives a choice of high-pass filters. The high-pass filter removes frequencies below the chosen frequency and allows high frequencies in the signal. These filters are useful for removing slowly moving baselines, such as motion or respiration artifacts, particularly in ECG (EKG) recordings.

Low-Pass Filtering. The Low Pass pop-up menu gives a choice of low-pass filters. The low-pass filter removes frequencies above the chosen frequency and allows low frequencies in the signal. These filters are useful for removing high-frequency signals, such as noise, and to prevent aliasing in the recorded signal.

Mains Filter. Select or deselect the Mains filter checkbox to turn the mains filter on and off. The mains filter allows you to remove interference related to the mains frequency (both fundamental and harmonic frequencies). This is an adaptive filter. It adjusts to filter the interference by tracking the input signal for a second. Because of this, in general, using the mains filter is a good option but it does have some limitations. More details on the mains filter can be found in the LabChart Help Center.

Inverting the Signal

The Invert checkbox provides a simple way to change the polarity of the recorded signal without having to swap the connections to the recording electrodes.

DC Restore

The DC Restore or Zero button is available in the Bio Amp dialog. It provides a quick way to bring the recording trace back to zero. It is for use with very low high-pass filter settings, since the trace can take a long time to zero automatically at those settings. In later versions of LabChart the **DC Restore All** command from the Setup menu performs this operation on all Bio Amp channels.

Using the Bio Amp

Standard electrophysiology texts describe how to record various types of biological signal. ADInstruments also produces materials describing specific uses of Bio Amps, such as Animal Physiology and Human Physiology LabChart experiments, that can be downloaded from the ADInstruments website or obtained from your ADInstruments representative.

Some Suitable Uses

The Bio Amp can measure a wide variety of biological signal sources. Some of these measurements include:

ECG. Electrocardiogram (also referred to as EKG); a recording of surface potentials due to electrical currents associated with the heartbeat.

EEG. Electroencephalogram; a recording of the electrical activity of the brain. Scalp electrodes record potential waves (10–100 μ V) representing the summed activity of cortical neurons.

EMG. Electromyography; a recording of the electrical activity of a muscle, using surface electrodes. The recorded activity may be a voluntary contraction, or evoked by motor nerve stimulation.

EOG. Electro-oculogram; a recording of the potential difference between the front and back of the eyeball, as projected on to the face. Only supported on Octal and Quad Bio Amps (not supported on Single and Dual Bio Amps)

ERG. Electroretinogram; a recording of the electrical signals produced in the retina by a light stimulus. Bilateral measurements require a Dual Bio Amp, or two Bio Amps.

Cortical Evoked Potentials. Averaged recordings of the electrical activity of the brain when subject to stimulation: visual evoked response, auditory evoked response, and somatosensory response. These should be done with signal averaging, using Scope view.

SNAP. Sensory nerve action potential; a recording of evoked responses in stimulated nerves. This is usually done with signal averaging, using Scope view.

Some Unsuitable Uses

The Single/Dual Bio Amp and Quad/Octal Bio Amp are not recommended for work requiring high-impedance electrodes or using a high bandwidth. Some of the tasks for which they are not really suitable include:

- Intracellular micropipette recordings. Recordings from a very fine, electrolyte-filled tube inserted into a nerve or muscle cell. These require an electrometer amplifier.
- Needle electromyography; the intramuscular recording of the electrical activity of a muscle, which requires low input capacitance and a driven guard.

For these a headstage is required. The bandwidth is limited to approximately 5 kHz.

Recording Technique

Several problems can arise when using the Bio Amp to record signals. These are basic problems of technique, and should be addressed before setting up. It is important to understand the types of problems that can occur, how they manifest, and what can be done to remove them or to minimize their effect. Potential problem areas include aliasing, frequency distortion, saturation, ground loops, electrode contact, motion artifacts, electromagnetic fields, and data display.

There is a good introduction to data acquisition provided in the documentation for LabChart. It is highly recommended reading for anyone recording biological signals using the front-end. Apart from the general areas covered in that material, two things particularly affect the kind of measurements made with the Bio Amp, and can cause ‘artifacts’ (spurious readings) in the recorded waveform: electrode contact and motion effects.

Bio Amp Operation

The Bio Amp and other ADInstruments front-ends have been designed to integrate fully into the PowerLab system. The Bio Amp is essentially an extension of the PowerLab’s input amplifiers. The amplification and ranges offered in LabChart result from the combination of both pieces of hardware.

The Bio Amp provides:

- full electrical isolation from power-line (mains) circuitry to guarantee subject safety
- a low-noise, high-gain differential amplifier specifically designed for biological signal measurements
- software-controlled low-pass, high-pass and adaptive mains filters to remove unwanted signal frequencies for particular uses
- audio output for use with EMG or EEG signals.

Right-leg Drive

All Bio Amps include a feature known as Right-leg drive or a Driven Right Leg circuit. This circuit helps to reduce electrical interference from the mains picked up by the human body or via patient leads. In order to use this feature, the patient connections of channel 1 must be connected to the subject and the Right-Leg drive checkbox must be ticked. If the Right-leg drive checkbox is unticked, the green common terminal is connected to isolated ground inside the Bio Amp. The Right-leg drive circuit uses the green Common terminal and a patient lead must be connected to this terminal in order to use this feature.

The Bio Amp Cable Input

The Bio Amp cable, supplied with the Bio Amp and Dual Bio Amp, plugs into the six-pin input socket on the front panel of the Bio Amp. A notch in the plug ensures that polarity is correct. Only the supplied Bio Amp cable and lead wires should be used as other cables may not meet safety requirements.

Note: The pin arrangements for the Bio Amp and the Dual Bio Amp are different; the Bio Amp cables are not interchangeable. (The Bio Amp is supplied with a 3-lead Bio Amp cable, whereas the Dual Bio Amp is supplied with a 5-lead Bio Amp cable.)

Troubleshooting

This section describes the common problems that can occur when using the Bio Amp with your PowerLab recording unit. It covers how these problems are caused, and what you can do to alleviate them. If the solutions here do not work, refer to earlier chapters, the LabChart Help Center, and your PowerLab Owner's guide for possible remedies. If none of the solutions here or elsewhere are of help, then consult your ADInstruments representative.

Most of the problems that users encounter are connection problems, and can usually be fixed by checking connections and starting up the hardware and software again. Very rarely will there be an actual problem with the front-end or the PowerLab.

Problems and Solutions

The status indicator fails to light when the software is started, or the Bio Amp... dialog(s) do not appear where they should

The I²C cable or one or more BNC-to-BNC cables from the Bio Amp to the PowerLab are not connected, have been connected incorrectly (to the wrong input, for instance), or are loose.

- Turn everything off. Check to see that all cables are firmly seated and screwed in. BNC cables from the Bio Amp must be connected to a positive input on the PowerLab. Make sure the input is the same channel from which you expect to use the front-end in the software. Start up again to see if this has fixed the problem.

You are using an early version of LabChart.

- Upgrade to the latest version of the software. Contact your ADInstruments representative for information.

The BNC or I²C cable is faulty.

- Replace the cable and try again. Immediately label all cables proved faulty so that you don't use them again by accident.

The Bio Amp is faulty

- This is the least likely event. If the front-end will not work properly after the previous measures, then try using it on another PowerLab. If the same problems recur with a second PowerLab, the front-end may be faulty. Contact your ADInstruments representative to arrange for repairs.

On starting up the software, an alert indicates that there is a problem with the front-end or driver

The correct Bio Amp driver is not installed on your computer.

- Reinstall the software.

You are using an early version of LabChart.

- Upgrade to the latest version of the software. Contact your ADInstruments representative for information.

The BNC or I²C cable is faulty.

- Replace the cable and try again. Immediately label all cables proved faulty so that you don't use them again by accident.

The Bio Amp is faulty.

- This is the least likely event. If the front-end will not work properly after the previous measures, then try using it on another PowerLab. If the same problems recur with a second PowerLab, the front-end may be faulty. Contact your ADInstruments representative to arrange for repairs.

The trace will not zero properly when adjusting high-pass filtering

The Bio Amp is receiving signals at a level that has saturated the input amplifier, resulting in a large offset. This is normally due to poor contact between the electrodes and the subject.

- Check the connections for proper contact and try again.

If you are using the Dual, Quad or Octal Bio Amps, with the 0.02 Hz high-pass filter, the trace can take a long time to zero by itself.

- Click the DC Restore button in the Bio Amp dialog to bring the trace back to zero quickly.

The signal appears to display a constant amplitude oscillation

Frequency interference from power lines can become superimposed on the biological signal being measured.

- You can use the mains filter in the Bio Amp to remove excessive line voltage frequency interference (use the checkbox in the Bio Amp dialog).

If you are using cables and leads that were not supplied with your Bio Amp, they may be unshielded or of low quality.

- Check to make sure that you are using high-quality shielded cables and high-quality leads. Only the supplied Bio Amp cable and leads should really be used.

The signal is noisy at lower ranges

This is probably the amplified noise from the electrodes, not a fault as such. There is, in addition, noise that cannot be avoided by any amplifier (called 'thermal' or 'Johnson' noise).

- Set the low-pass filter to remove the noise. (But be careful, since important components of the signal could also be attenuated.)

This could be due an electrically noisy environment, particularly if there is some equipment that produces a radio frequency that interacts with the Bio Amp modulator, giving a heterodyne effect.

-
- Turn off pieces of unnecessary equipment to try and isolate the cause, then either leave the equipment off, or, if possible, move the subject or equipment outside the area of any interfering field.
 - At the lowest ranges you may have to shield, shorten, or even replace the Bio Amp leads, since they will tend to act as radio receptors.

The signal appears to be unusual, very weak, clipped, or distorted in some way

This may be a problem of technique: the sampling rate, range, or filter settings may be inappropriate for the signal you are recording.

- Make sure the settings are appropriate for the expected signal.

Specifications

Single Bio Amp

Input

Connection type:	Three shrouded 1.5 mm male pin sockets
Input configuration:	1 isolated differential channel with isolated ground reference or right leg drive
Input impedance:	>1 G Ω differential, <100 pF (no cable) or 500 pF (supplied shielded Bio Amp cable and leads) to isolated ground

Physical Configuration

Dimensions (h \times w \times d):	55 mm \times 120 mm \times 260 mm (2.2" \times 4.7" \times 10.2")
Weight:	1.3 kg (2 lb 12 oz)
Power requirements:	~2 W

Dual Bio Amp

Input

Connection type:	Six-pin DIN/MS socket to fit 5-lead Bio Amp cable (Tronomed D-1540)
Configuration:	2 isolated differential channels with common isolated ground reference or right leg drive
Input impedance:	>1 G Ω differential, <100 pF (no cable) to isolated ground or 500 pF (supplied shielded Bio Amp cable and leads) to isolated ground

Physical Configuration

Dimensions (h \times w \times d):	55 mm \times 120 mm \times 260 mm (2.2" \times 4.7" \times 10.2")
Weight:	1.3 kg (2 lb 12 oz)
Power requirements:	~3 W

Quad/Octal Bio Amp

Input

Connection type:	9 \times 1.5 mm (Quad) or 17 \times 1.5 mm (Octal) pin shrouded male socket to suit single pin 1.5 mm
------------------	---

Configuration:	4 (Quad) or 8 (Octal) isolated differential channels with common isolated ground reference or right leg drive
Input impedance:	>1 G Ω differential, <100 pF (no cable) to isolated ground

Physical Configuration

Dimensions (h × w × d):	70mm x 240mm x 260mm (2.7" x 9.45" x 10.2")
Weight:	2.5 kg (5.5 lb)
Power requirements:	~4.5W (Quad), ~7W (Octal)

Common specifications for all Bio Amps

Input

Isolation:	4000 V _{rms} (50 Hz for 1 minute)
Input ranges:	$\pm 100 \mu\text{V}$ to $\pm 100 \text{ mV}$ full scale in 10 steps (combined PowerLab and Bio Amp) $\pm 100 \text{ mV}$ $\pm 50 \text{ mV}$ $\pm 20 \text{ mV}$ $\pm 10 \text{ mV}$ $\pm 5 \text{ mV}$ $\pm 2 \text{ mV}$ $\pm 1 \text{ mV}$ $\pm 500 \mu\text{V}$ $\pm 200 \mu\text{V}$ $\pm 100 \mu\text{V}$
Gain accuracy:	$\pm 1.5\%$ all ranges
Non-linearity:	< 0.2% within range
Noise at various bandwidths:	1 Hz to 10 kHz: < 1.2 μV_{rms} (< 12 μV p-p) 0.3 Hz to 1 kHz: < 0.5 μV_{rms} 0.1 Hz to 100 Hz: < 0.4 μV_{rms}
IMRR (isolation mode):	> 130 dB (to non-isolated earth, 50 Hz)
CMRR (common mode):	>60 dB typical (150K electrode impedance, 5K imbalance @50Hz and 60Hz) >100dB typical (Balanced electrode impedance @50Hz and 60Hz)

Input leakage current: $< 4 \mu A_{rms}$ @ 240 V, 50 Hz

DC tolerance: ± 370 mV

Baseline restoration: Automatic or manual

Filtering

Low-pass filtering: Frequencies software-selectable.

50, 100, 200, 500, 1000, 2000, 5000, 10000 Hz, Off

High-pass filtering: Frequencies software-selectable.

D.C, 0.0003, 0.001, 0.003, 0.01, 0.03, 0.1, 0.3, 1, 3, 10, 30, 100, 200, 300 Hz

Adaptive mains filter: 50 or 60 Hz frequency (Refer to Powerlab owner's guide for further information)

Output

Signal: ± 4.5 V maximum

Audio output: 3.5mm stereo output socket suitable for direct headphone or powered speaker connection. Output selectable from software.

Control Port

I²C port: Provides control and power. Interface communications rate of ~50 kbits/s.

Physical Configuration

Operating conditions: 5–35 °C, 0–90% humidity (non-condensing)

Transport/Storage conditions: 0–40 °C, 0–95% humidity (non-condensing)



NOTE:

PowerLabs provide power to the Bio Amp and no other power source is required.

Regulatory Information

Safety:	Complies with IEC 60601-1:2012 (tested by TUV Singapore)
EMC:	Complies with IEC 60601-1-2:2014 (tested by EMC Technologies, Sydney, Australia)
Equipment:	Use only with an ADInstruments 35 series PowerLab. The PowerLab must be connected to safety earth via the power supply cable to ensure electrical safety.
Operation:	Continuous

Body protection rating

(Applied parts only):



(BF)

Unsuitable uses:	Do not use in the presence of flammable anaesthetic - air mixtures. Avoid operating near high voltage, RF or strong magnetic fields that may cause interference.
Method of Disposal:	Forward to recycling centre or return to manufacturer.

This equipment is not intended to be modified or serviced by the user. No user serviceable parts inside. Refer servicing to authorised ADInstruments service centre. ADInstruments reserves the right to alter these specifications at any time.

Electromagnetic Compatibility

The Bio Amps (the devices) have been tested to comply with the requirements of IEC 60601-1-2. (The relevant basic EMC standards are IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-6, IEC 61000-4-8, IEC 61000-4-11 and CISPR 11).

Emissions

- The devices use RF energy for its internal function only. RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
- The devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Immunity

- Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
- Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
- Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Separation Distances

- The devices are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.
- Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance in the table below.
- Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.

Rated maximum output power of transmitter, P	Separation distance	
	150 kHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1.17\sqrt{P}$	$d = 2.33\sqrt{P}$
0.01 W	0.1 m	0.2 m
0.1 W	0.4 m	0.7 m
1 W	1.2 m	2.3 m
10 W	3.7 m	7.4 m
100 W	11.7 m	23.4 m



Chapter 4

Neuro Amp EX

WARNING:
Refer to Intended
Use statement
on page 1 of this
Owner's Guide
before use.

The FE285 Neuro Amp EX is one of a family of devices called front-ends, designed to extend the capabilities of the PowerLab system. It differs from other front-ends by including an additional piece of hardware, the headstage. The Neuro Amp EX allows the PowerLab system to record extracellular action potentials from single cells or groups of cells.

The Neuro Amp EX provides:

- full electrical isolation from power-line (mains) circuitry to guarantee subject safety.
- a low-noise, high-gain differential amplifier specifically designed for neurophysiological extracellular signal measurements.
- software-controlled low-pass, high-pass and mains filters to remove unwanted signal frequencies for particular uses.
- audio output to listen to neurological signals.



The Neuro Amp EX is fully isolated for human use, in addition to being suitable for making extracellular recordings from experimental animals.

The Neuro Amp

The Neuro Amp EX [FE285] is designed to allow the PowerLab system to perform isolated measurements of single-unit or multi-unit extracellular action potentials from single neurons or axons, from human or animal subjects. It consists of an electrically isolated differential input AC amplifier with an isolated ground connection.


It is possible to 'daisy-chain' multiple Neuro Amp EXs to record from separate subjects, using separate grounds. However, it is not recommended to use more than one Neuro Amp EX to record from a single subject or to use a Neuro Amp EX with an ADInstruments' Bio Amp to record from a single subject (this is discussed in more detail in the next section).

The Front Panel

The front panel of a Neuro Amp EX has an input socket and a status indicator light.

Figure 4-1

The front panel of the Neuro Amp EX



The Input Socket

The headstage output cable is coupled to the Neuro Amp EX front-end using a five-pin input socket on the front panel. The pin arrangement ensures that the polarity is correct. The input socket is physically and electrically isolated from the low-voltage mains-supply circuitry of the PowerLab and the input connections are isolated internally by isolation circuitry. The socket provides 7.5 V supply lines to the headstage, a protected earth and differential input lines.

The Status Indicator

The status indicator light of a Neuro Amp EX is located at the bottom left of the front panel.

When an ADInstruments application such as LabChart starts, the status indicator should flash briefly and then remain green, indicating that the program has found the front-end, checked and selected it, and is ready to use it. If the status indicator does not turn on and stay on when the application starts, it is most likely that the front-end is not connected properly.

WARNING:
When used in ambient temperatures of 38 degrees C and above, do not touch the Neuro Amp enclosure for more than a minute continuously.

The Back Panel

The back panel of the Neuro Amp EX provides all the sockets required to connect the front-end to the PowerLab and to other front-ends.




Figure 4-2

The back panel of the Neuro Amp EX

WARNING:

The Neuro Amp has only been assessed as compliant with IEC60601-1 safety standard when used with a 35 series PowerLab

I²C Input and Output Sockets

Two nine-pin sockets are used to communicate with the PowerLab (they are marked 'I²C Bus': a 'bus' is simply information-transmission circuitry such as cables and connectors). These sockets allow multiple front-ends to be used independently with one PowerLab.

Power and control signals to the front-ends come from the PowerLab. Many front-ends can be connected to the system, in series, output to input, providing there is the same number of channel inputs available on the PowerLab (this is discussed in more detail in the next chapter).

Analog Out Socket

The BNC socket labeled Analog Out on the back panel provides the signal output and a connection is made to an analog input socket on the front of the PowerLab with a BNC-to-BNC cable (supplied).

If you are connecting to an earlier model PowerLab with differential inputs, only connect the analog output to the positive analog input of the PowerLab. ADInstruments applications do not find a front-end on start up if the negative input is used.

Audio Out Socket

The Neuro Amp EX has an audio monitor output on the back panel that can be used with a wide range of headphones or externally powered speakers. The 3.5 mm socket provides mono sound that may be of use when monitoring nerve firings, to control the placement of electrodes for instance.

WARNING:

As with all audio devices, there is a risk of temporary hearing damage if used carelessly. It is recommended headphones are attached to the Neuro Amp or Bio Amp audio socket prior to attaching the headphones to ears..

Using More Than One Neuro Amp EX

Although you can record with more than one Neuro Amp EX front-end using a common Ground/reference, such as when recording from a single subject, it is not recommended. If using two or more Neuro Amp EX front-ends, some interaction can take place between them, causing up to 10 μV of induced low-frequency signal, due to slight differences in frequency between the isolated power supplies. This should not affect most biopotential measurements, especially if using signal averaging. Narrowing the signal bandpass may also reduce the effect.

Software Requirements

The Neuro Amp EX requires the following versions of ADInstruments software applications:

- LabChart version 8.1.11 and above

Note: the Neuro Amp EX may not operate correctly with earlier versions of these applications. If you have queries regarding hardware and software requirements of the Bridge Amps, please contact your local ADInstruments representative.

Using LabChart

When a Neuro Amp EX is properly connected to the PowerLab, the **Input Amplifier...** menu command is replaced with **Neuro Amp EX...** for the input to which the Neuro Amp EX is connected. If the application fails to find a front-end connected, the normal text remains. If the expected text is not displayed, quit the application, turn the PowerLab off and check the BNC and I²C connections. Then restart the PowerLab and relaunch the application.


The documentation for LabChart does not cover front-end-specific features. These features are described in detail here for LabChart. Generally, dialogs for these applications are similar. The Neuro Amp EX dialog for LabChart for Macintosh and LabChart for Windows are very similar and are described here together.

Neuro Amp EX dialog

The Neuro Amp EX dialog (Figure 9–5 and Figure 9–6) allows software control of the combined input amplifiers and filters in the PowerLab and Neuro Amp EX. The signal present at a channel's input is displayed in the preview area so that you can see the effects of changes in settings. After changing settings in the dialog, click the **OK** button to apply them.

The Neuro Amp EX dialog appears when you choose **Neuro Amp EX...** from a Channel Function pop-up menu. To set up many channels quickly, open the **Setup > Channel Settings...** dialog. Here you can view all the channels that are turned on, and you can turn off any unnecessary channels. Clicking on **Neuro Amp EX...** in the Input Settings column of the Channel Settings dialog will also open the Neuro Amp EX dialog.

Figure 4-3
The Neuro Amp
EX dialog for
Windows



Signal Display

The input signal is displayed so you can see the effect of changing the settings – no data is recorded while setting things up. The average signal value is shown above the display area.

You can shift and stretch the vertical Amplitude axis by clicking and dragging it in various ways to make the best use of the available display area. It functions the same as the Amplitude axis of the Chart Window; the controls are identical and any change is also applied to the Chart Window.

Setting the Range

The Range pop-up menu lets you select the input range or sensitivity of the channel – the combined range of the PowerLab and the Neuro Amp EX. Changing the range in the Neuro Amp EX dialog is equivalent to changing it in the LabChart window. The default setting is 1 mV and the ranges go down to 20 µV in 6 steps.

Filtering

The Neuro Amp EX has low-pass, high-pass and mains filter circuitry that can be adjusted to suit the application. The mains filter removes mains frequency related interference. In general, it is better to prevent interference at the source than to filter it. The high-pass and low-pass filters provide bandwidth limiting of low-frequency and high-frequency

signals, respectively. Note that the settings for one filter type may restrict the possible settings for the other.

High-Pass Filtering. The High Pass pop-up menu gives the choice of 100 Hz and 300 Hz high-pass filters. The high-pass filter removes frequencies below the selected cut-off frequency and allows high frequencies in the signal to pass. Note: this amplifier is not suitable for DC recording.

Low-Pass Filtering. The Low Pass pop-up menu gives the choice of 1 kHz, 2 kHz and 5 kHz low-pass filters. The low-pass filter removes high-frequency components above the selected cut-off frequency and allows lower frequencies in the signal to pass. These filters are useful for removing high-frequency signals, such as noise, and to prevent aliasing in the recorded signal (see page 24).

Mains Filter. Select or deselect the Mains filter checkbox to turn the mains filter on and off. The mains filter allows you to remove interference related to the mains frequency (both fundamental and harmonic frequencies). This is an adaptive filter. It adjusts to filter the interference by tracking the input signal for a second, creating a template of the interfering frequencies and then subtracting this template from the input signal. Because of this, in general, using the mains filter is a good option. However, the mains filter does have some limitations, such as not being useful for very short recordings of less than one second. More details on the mains filter can be found in the **LabChart Help Center**.

Inverting the Signal

The Invert checkbox provides a simple way to change the polarity of the recorded signal without having to swap the connections to the recording electrodes.

Types of Measurement

The low noise and high gain of the Neuro Amp EX makes it suitable for all recordings that require a high bandpass (300 Hz – 5 kHz) and a high signal-to-noise ratio. Such recordings may include extracellular recordings from single cells or groups of cells, or from nerve fibres (split-nerve preparation or human microneurography). The headstage provides a gain of 100x and the cable shielding is directly connected to the casing, limiting the need for additional shielding at the input terminals.


Note: The Neuro Amp EX is not suitable for intracellular recording with glass microelectrodes.

The Neuro Amp EX Headstage

The headstage can be mounted in a micromanipulator for animal use or taped directly to the skin of a human subject. The casing of the headstage is connected to the shielding of the cable, thereby providing shielding up to the input terminals. Three female gold plated miniature sockets provide inputs on the headstage. The positive and negative


Figure 4-4
Input sockets for
the Neuro Amp EX
headstage

terminals are adjacent, while the ground terminal is further away. The positive terminal is distinguished by a slightly recessed surround.



Connections are made to the Neuro Amp EX front-end using the supplied Neuro Amp EX headstage and connectors (microelectrodes are not supplied). The headstage cable plugs into the five-pin input socket on the front panel: the pin arrangement ensures that polarity is correct. Only use the supplied Neuro Amp EX headstage as others may not meet safety requirements.

Figure 4-5
Headstage for
the Neuro Amp
EX, as supplied
(top), and with
male connectors
attached (below)



Six male gold plated miniature connectors or pins are supplied with the headstage. These can be crimped or soldered to standard metal microelectrodes, such as tungsten, stainless steel or platinum, to allow you to construct electrodes to suit your requirements.

Recording Technique

When conducting research using the Neuro Amp EX, it is important to ensure that you use correct preparation and recording techniques. It is important to understand the types of problems that can occur, how they manifest, and what can be done to remove them or to minimize their effect. Potential problems may arise due to aliasing, frequency distortion, saturation, ground loops, electromagnetic fields, electrode contact and motion artifacts. Many of these are discussed in “Preventing Problems” on page 14. See the section below for discussion of preventing artifacts due to poor electrode contact and subject movement.

Setting up to Record From a Subject

This section provides a step-by-step guide to making a recording from a subject. After making the subject comfortable, a suitably qualified operator can use the following steps to attach the headstage. The operator should follow their institute's guidelines for making neurological recordings and use precautions to minimize the risk of infection:

1. Attach the headstage to the Neuro Amp EX.
2. Securely attach the headstage to the subject by taping or bandaging it to the limb or trunk.
3. Insert all the pins into the headstage. Attach the Ground/reference electrode first, using Figure 9–7 as a guide.
4. Connect the active microelectrode (+) and the reference electrode (-) to the positive and negative inputs of the headstage, respectively.
5. The Audio Output from the Neuro Amp EX can be connected to external speakers, if required.
6. To help locate the nerve site, nerve stimulating techniques compatible with your institute's guidelines may be used. Stimulation may be made using the ADInstruments Stimulus Isolator, or another Isolated Stimulator.
7. Changes in signal output in response to nerve stimulation can be observed using the Neuro Amp EX dialog or by recording data using Scope View in LabChart.
8. If no longer required, move stimulation leads away from the headstage before recording. Leaving them attached or close to the headstage will introduce noise.


For microneurography studies, the subject is normally passive during the procedure. The headstage *must* be connected to the Neuro Amp EX before and during the procedure. During the procedure, an operator:

- securely attaches the headstage to the subject using tape or a bandage (Figure 9–9);
- chooses a zero/reference point on the subject by attaching the Ground/Reference electrode to the subject (Figure 9–10);
- and, finally, connects the active and reference microelectrodes to the subject (Figure 9–11).

The headstage must be securely attached (Figure 9–9) to the limb of the subject under investigation. This ensures that the headstage is in contact with the subject and helps to minimize artifacts from movement.

Figure 4–6

The headstage should be securely attached to the subject.



The Ground/Reference lead can be attached to a surface electrode using a press stud or snap fitting (Figure 9–10).

Figure 4–7


The Ground/Reference lead wire attached with a press stud



Microelectrode needles should be positioned using examination gloves to minimize the risk of infection. The active microelectrode is inserted into the subject's skin overlying the nerve, while the reference electrode is inserted under adjacent skin. Under no circumstances should the headstage be disconnected from the Neuro Amp EX while the electrodes are still attached to the subject. Not having the headstage attached to the Neuro Amp EX before and during the procedure may result in small static electric shocks being delivered to the subject through the microelectrodes.

Figure 4–8

The active microelectrode is inserted over the common peroneal nerve at the fibular head. The reference electrode is inserted under the skin. The ground (Ag/AgCl) electrode is attached to the surface of the skin with the green lead.



Motion Effects

A common source of artifacts when recording neurological signals is motion of the subject or equipment. For example, muscular activity generates its own electrical signals, which may interfere with the neuronal signal, depending on the location of the electrodes.

If the subject is going to move during recording, then special care needs to be taken when positioning electrodes and securing leads.

Electrode Contact

Occasionally during measurement of a neurological signal, one of the lead wires connecting the source to the front-end may become disconnected, or electrode contact may become poor. If this should happen, relatively high voltages (potentials) can be induced in the open wire, due to electric fields caused by the power line or other sources close to the front-end or to the subject. These induced potentials result in a constant amplitude disturbance in the recorded signal at the mains frequency, causing masking or loss of the desired signal. If the problem is a recurring one, one of the leads may be faulty. Check connections and replace faulty leads, if necessary.


Technical Aspects

The Neuro Amp EX and other ADInstruments front-ends have been designed to integrate fully into the PowerLab system. The Neuro Amp EX is essentially an extension of the PowerLab's input amplifiers. The amplification and ranges offered in LabChart result from the combination of both pieces of hardware.

The PowerLab provides control and low-voltage power to front-ends through a special expansion connector called the I²C bus. Front-ends are also connected to the analog inputs of the PowerLab via a BNC cable, through which they send the amplified and filtered signals. The overall operation of the Neuro Amp EX can be better understood by referring to Figure 9–13 on page 51.

The headstage is a differential in/differential out amplifier with a fixed gain of $\times 100$ and a 10 Hz high-pass filter. The amplifier drives a differential signal into the Neuro Amp EX. Power is supplied through the connector. A reference pin is a current-limited isolated ground return to protect against excessive auxiliary current. The overall operation of the Neuro Amp EX headstage can be better understood by referring to Figure 9–12.

Figure 4–9
Block diagram
of the headstage
amplifier




An audio signal output, capable of driving headphones or powered speakers, is provided by tapping off and buffering the output stage.

The control for the various filters and gain stages in the Neuro Amp EX is provided by on-board microprocessors, which also communicate with the PowerLab over the I²C bus.

The Neuro Amp EX Input

The Neuro Amp EX headstage cable plugs into the five-pin input socket on the front panel of the Neuro Amp EX (Figure 9–14). The pin arrangement ensures that polarity is correct. Only the supplied Neuro Amp EX headstage should be used. Using other headstages invalidates the warranty and may not be as safe. ADInstruments are not responsible for damages incurred if using third party headstages.

Figure 4–10
Pin assignments
for the Neuro
Amp EX input
connector



The Neuro Amp EX has one connector: the socket provides two pins for a differential input signal, one pin for the ground/reference, and two pins for power supply to the headstage. The entire connector is physically and electrically isolated to ensure subject safety.

Applied parts of
the Neuro Amp



Supplied Accessories (Applied parts)

- MLT185 Neuro Amp Headstage
- WPI 220-883-S02 miniature sockets.

Troubleshooting

This appendix describes most of the common problems that can occur when using the Neuro Amp EX with your PowerLab recording unit. It covers how these problems are caused, and what you can do to alleviate them. If the solutions here do not work, earlier chapters, the LabChart Help Center, and the guide to your PowerLab may contain possible solutions. If none of the solutions here or elsewhere are of help, then consult your ADInstruments representative.

Most of the problems you may encounter are connection problems, and can usually be fixed by checking connections and starting up the hardware and software again. Very rarely will there be an actual problem with the front-end or the PowerLab.

Problems and Solutions

The status indicator fails to illuminate when the software is started, or the front-end commands do not appear where they should

The I²C cable or the BNC cable from the front-end to the PowerLab are not connected, have been connected incorrectly (to the wrong input, for instance) or are loose.

- Turn everything off. Check to see that all cables are firmly seated and screwed in. The BNC cable from the Neuro Amp EX must be connected to a positive input on the PowerLab. Make sure the input is the same channel from which you expect to use the front-end in the software. Start up again to see if this has fixed the problem.

You are using an early version of LabChart.

- Upgrade to the latest version of the software. Contact your ADInstruments representative for information.

The BNC or I²C cable is faulty.

- Replace the cable and try again. Immediately label all cables proved faulty so that you don't use them again by accident.

The Neuro Amp EX is faulty.

- This is the least likely event. If the front-end will not work properly after the previous measures, try using it on another PowerLab. If the same problems recur with a second PowerLab, the front-end may be faulty. Contact your ADInstruments representative to arrange for repairs.

On starting up the software, an alert indicates that there is a problem with the front-end or driver

The correct driver is not installed on your computer.

- Reinstall the software.

You are using an early version of LabChart.

- Upgrade to the latest version of the software. Contact your ADInstruments representative for information.

The BNC or I²C cable is faulty.

- Replace the cable and try again. Immediately label all cables that proved faulty so that you don't use them again by accident.

The Neuro Amp EX is faulty.

- This is the least likely event. If the front-end will not work properly after the previous measures, try using it on another PowerLab. If the same problems recur with a second PowerLab, the front-end may be faulty. Contact your ADInstruments representative to arrange for repairs.

The trace will not zero properly when adjusting high-pass filtering

The Neuro Amp EX is receiving signals at a level that has saturated the input amplifier, resulting in a large offset. This is normally due to poor contact between the electrodes and the subject.

- Check the connections for proper contact and try again.

The signal appears to display a constant amplitude oscillation

Frequency interference from power lines can become superimposed on the biological signal being measured.

- You can use the mains filter in the Neuro Amp EX to remove excessive line voltage frequency interference (use the checkbox in the Neuro Amp EX dialog).

The leads you are using may be unshielded or of low quality.

- Check to make sure that you are using high-quality leads. Only the supplied Neuro Amp EX headstage cable should be used.

The signal is noisy at lower ranges

This is probably the amplified noise from the electrodes, not a fault as such. There is, in addition, noise that cannot be avoided by any amplifier (called 'thermal' or 'Johnson' noise).

- Set the low-pass filter to remove the noise. (But be careful, since important components of the signal could also be attenuated.)

This could be due an electrically noisy environment, particularly if there is some equipment that produces a radio frequency that interacts with the Neuro Amp EX modulator, giving a heterodyne effect.

- Turn off pieces of unnecessary equipment to try and isolate the cause, then either leave the equipment off or, if possible, move the subject or equipment outside the area of any interfering field.
- At the lowest ranges you may have to shield, shorten, or even replace the leads, since they will tend to act as radio receptors.

The signal appears to be unusual, very weak, clipped or distorted in some way

This may be a problem of technique: the sampling rate, range or filter settings may be inappropriate for the signal you are recording.

- Make sure the settings are appropriate for the expected signal.

Specifications

Neuro Amp EX Front-end [FE285]

Input

Connection type:	Five-pin Redel connector	
Configuration:	One isolated differential channel with isolated ground reference	
Input impedance:	100 M Ω	
Safety:	Approved to IEC 60601-1 Standard (BF rating)	
Isolation:	4000 V _{rms} (50 Hz for 1 minute)	
Input ranges:	$\pm 20 \mu\text{V}$ to $\pm 1 \text{ mV}$ full scale in 6 steps (combined PowerLab, Neuro Amp EX front-end and headstage)	
Gain:	Range	Resolution
	1 mV	500 nV
	500 μV	250 nV
	200 μV	100 nV
	100 μV	50 nV
	50 μV	25 nV
	20 μV	10 nV

Filtering

Low-pass filtering:	Software-selectable: 1 kHz, 2 kHz, 5 kHz.
High-pass filtering:	First-order filter. Software-selectable: 100 Hz, 300 Hz.

Output

Signal:	$\pm 4.5\text{V}$ Maximum
Audio:	Stereo jack with mono output; $\pm 200 \text{ mV}$ full scale, current limited to $\pm 5 \text{ mA}$. Suitable for headphones or powered speakers.

Control Port

I ² C port:	Provides control and power.
------------------------	-----------------------------



NOTE:
PowerLabs provide power to the Neuro Amp and no other power source is required.

Physical Configuration

Dimensions (h × w × d):	55 mm × 240 mm × 260 mm (2.2" × 4.7" × 10.2")
Weight:	1.3 kg (2 lb 13 oz)
Power requirements:	~2 W
Operating conditions:	5–35 °C, 0–90% humidity (non-condensing)
Transport/Storage conditions:	0–40 °C, 0–95% humidity (non-condensing)

Regulatory Information

Safety:	Complies with IEC 60601-1:2012 (tested by TUV Singapore)
EMC:	Complies with IEC 60601-1-2:2014 (tested by EMC Technologies, Sydney, Australia)
Equipment:	Use only with an ADInstruments 35 series PowerLab. The PowerLab must be connected to safety earth via the power supply cable to ensure electrical safety.
Operation:	Continuous

Body protection rating

(Applied parts only):



(BF)

Unsuitable uses:	Do not use in the presence of flammable anaesthetic - air mixtures. Avoid operating near high voltage, RF or strong magnetic fields that may cause interference.
Method of Disposal:	Forward to recycling centre or return to manufacturer.

Neuro Amp EX Headstage [MLT185]

Connection type:	Gold plated miniature sockets WPI 220-883-S02. To suit Male 220-P02
Configuration:	Differential
Input impedance:	100 MΩ // ~100 pF
Gain:	×100
Input leakage current:	~1 pA
CMRR:	100 dB typical @ 50/60 Hz
Noise (rti):	< 2 μV _{rms} , < 14 μV _{pp} (100 Hz – 5 kHz)
Power:	From the Neuro Amp EX front-end

Electrode type:	Suitable for metal microelectrodes (300 kΩ typical)
Dimensions:	8 mm OD, 100 mm long
Material:	Stainless steel
Cable length:	2.8 m

This equipment is not intended to be modified or serviced by the user. No user serviceable parts inside. Refer servicing to authorised ADInstruments service centre. ADInstruments reserves the right to alter these specifications at any time.

Electromagnetic Compatibility

The FE285 Neuro Amp EX (the device) has been tested to comply with the requirements of IEC 60601-1-2, (IEC 61000-3-2, IEC 61000-3-3, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-8, IEC 61000-4-11 and CISPR 11).

Emissions

- The device uses RF energy for its internal function only. RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
- The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Immunity

- Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
- Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
- Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Separation Distances

- The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.
- Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance in the table below.
- Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.

Rated maximum output power of transmitter, P	Separation distance	
	150 kHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.17\sqrt{P}$	$d = 2.33\sqrt{P}$
0.01 W	0.1 m	0.2 m
0.1 W	0.4 m	0.7 m
1 W	1.2 m	2.3 m
10 W	3.7 m	7.4 m
100 W	11.7 m	23.4 m



Chapter 1

Warranty

Product Purchase and License Agreement

This Agreement is between ADInstruments NZ Ltd [‘ADI’] and the purchaser [‘the Purchaser’] of any ADI product or solution — software, hardware or both — and covers all obligations and liabilities on the part of ADI, the Purchaser, and other users of the product. The Purchaser (or any user) accepts the terms of this Agreement by using the product or solution. Any changes to this Agreement must be recorded in writing and have ADI’s and the Purchaser’s consent.

Responsibilities

The Purchaser and any others using any ADI product or solution agree to use it in a sensible manner for purposes for which it is suited, and agree to take responsibility for their actions and the results of their actions. If problems arise with an ADI product, ADI will make all reasonable efforts to rectify them. This service may incur a charge, depending on the nature of the problems, and is subject to the other conditions in this Agreement. ADI does not separately warrant the performance of products, equipment or software manufactured by third parties which may be provided to Purchaser as part of an overall solution. However, as further noted below, ADI will pass through to Purchaser all applicable third party warranties to the extent it has the right to do so.

ADI Product Hardware Warranty

ADI warrants that PowerLab Data Acquisition Units (PL prefix)¹ and Front-ends (FEprefix)² shall be free from defects in materials and workmanship for five (5) years from the date of purchase. Other PowerLab Data Acquisition Units³, Front-ends⁴ and Pods⁵ shall be free of defects in material and workmanship for three (3) years from their date of purchase. ADI also warrants that ADI Specialized Data Recorders⁶ and Instruments⁷ shall be free of defects in material and workmanship for one (1) year from their date of purchase. If there is such a defect, as Purchaser’s sole remedy hereunder, ADI will repair or replace the equipment as appropriate, and the duration of the warranty shall be extended by the length of time needed for repair or replacement.

To obtain service under this warranty, the Purchaser must notify the nearest ADI office, or Authorized Representative, of the defect before the warranty expires. The ADI or Representative office will advise the Purchaser of the nearest service center address to which the Purchaser must ship the defective product at his or her own expense. The product should be packed safely, preferably in its original packaging. ADI will pay return shipping costs.

Hardware Warranty Limitations

This warranty applies only to the ADI hardware specified in this document and used under normal operating conditions and within specification. Consumables, electrodes and accessories are not covered by this warranty. Third party equipment may be covered by the third party manufacturer's warranty. To the extent that ADI has the right to pass through any third party manufacturer warranties to Purchaser it will do so to the extent it is able to do so. Copies of applicable third party manufacturer warranties, to the extent they exist, are available upon request. The warranty provided hereunder does not cover hardware modified in any way, subjected to unusual physical, electrical or environmental stress, used with incorrectly wired or substandard connectors or cables, or with the original identification marks altered. Tampering with or breaking of the Warranty Seal will also void the warranty.

Product Types & Warranty Term

ADI manufactured products covered by a five (5) year warranty

¹ Data Acquisition Units: PowerLab 35 series with PL prefix

² Front-ends: ADI Front-end Signal Conditioners with FE prefix.

ADI manufactured products covered by three (3) year warranty

³ Data Acquisition Units: PowerLab 26 series with ML prefix

⁴ Front-ends: ADI Front-end Signal Conditioners with ML prefix.

⁵ Pods: The entire range of ADI Pod Signal Conditioners.

ADI manufactured products covered by one (1) year warranty

⁶ Specialized Data Recorders: Metabolic Systems (e.g., ML240 PowerLab/8M Metabolic System)

⁷ Instruments: Blood FlowMeter, Gas Analyzers, NIBP System (excluding transducers), STH Pump Controller.

Third Party Products (Including Transducers)

Products not manufactured by ADI are covered by the manufacturer's warranty.

Accessories and Consumables

Accessories and Consumables are not covered by any type of warranty.

General Limitations

ADI products are produced to high standards, and should perform as described in the supplied documentation. There is a limited hardware warranty, and technical support is provided for all ADI products. Nevertheless, since ADI products could be affected by external factors (for instance, the computer system on which they run and other hardware and/or software provided by third parties), absolute performance and reliability of products and the overall solution cannot be guaranteed. No warranty, either expressed or implied or statutory, other than that expressly contained in this Agreement, is made in respect to ADI products or software, third party products or software, the overall solution or otherwise. The Purchaser therefore assumes all risks as to the performance and reliability of the products, the software, the solution and the results gained using them. ADI neither assumes or authorizes any person to assume on its behalf any liability in connection with the sale, installation, service or use of its products. ADI shall not be held responsible for special, consequential or punitive damages of any kind arising out of

sale, installation service or use of its products.

EXCEPT FOR THE EXPRESS WARRANTY SET FORTH HEREIN, THE SOLUTION AS WELL AS ALL EQUIPMENT AND SOFTWARE PROVIDED HEREUNDER ARE PROVIDED “AS IS” AND ADI MAKES NO WARRANTY. AS TO ITS USE OR PERFORMANCE. EXCEPT FOR ANY WARRANTY, CONDITION, REPRESENTATION OR TERM THE EXTENT TO WHICH CANNOT BE EXCLUDED OR LIMITED BY APPLICABLE LAW, ADI AND ITS SUPPLIERS MAKE NO WARRANTY, CONDITION, REPRESENTATION, OR TERM (EXPRESS OR IMPLIED, WHETHER BY STATUTE, COMMON LAW, CUSTOM, USAGE OR OTHERWISE) AS TO ANY MATTER INCLUDING, WITHOUT LIMITATION, NON INFRINGEMENT OF THIRD PARTY RIGHTS, MERCHANTABILITY, INTEGRATION, OR FITNESS FOR A PARTICULAR PURPOSE. YOU ASSUME RESPONSIBILITY FOR SELECTING THE SOLUTION TO ACHIEVE YOUR INTENDED RESULTS, AND FOR THE INSTALLATION OF, USE OF, AND RESULTS OBTAINED FROM THE EQUIPMENT AND SOFTWARE. WITHOUT LIMITING THE FOREGOING PROVISIONS, ADI MAKES NO WARRANTY THAT THE EQUIPMENT OR SOFTWARE WILL BE ERROR-FREE OR FREE FROM INTERRUPTIONS OR OTHER FAILURES OR THAT THE SOFTWARE OR EQUIPMENT WILL MEET YOUR REQUIREMENTS. UNDER NO CIRCUMSTANCES AND UNDER NO LEGAL THEORY, WHETHER IN TORT, CONTRACT, OR OTHERWISE, SHALL ADI OR ITS SUPPLIERS BE LIABLE TO PURCHASER OR TO ANY OTHER PERSON FOR LOSS OF PROFITS, LOSS OF GOODWILL, OR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, OR DAMAGES FOR GROSS NEGLIGENCE OF ANY CHARACTER INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF GOODWILL, WORK STOPPAGE, COMPUTER FAILURE OR MALFUNCTION, OR FOR ANY OTHER DAMAGE OR LOSS. IN NO EVENT SHALL ADI OR ITS SUPPLIERS BE LIABLE FOR ANY DAMAGES IN EXCESS OF THE PRICE PAID FOR THE EQUIPMENT AND SOFTWARE, EVEN IF ADI, OR ITS AUTHORIZED PARTNERS OR SUPPLIERS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

ADI is acting on behalf of its suppliers for the purpose of disclaiming, excluding and/or limiting obligations, warranties and liability as provided in this agreement, but in no other respects and for no other purpose. The foregoing provisions shall be enforceable to the maximum extent permitted by applicable law.

Controlling Law and Severability

This license shall be governed by the laws of the territory into which the software is sold, or if sold into the United States of America, by the laws of the State of California.

Technical Support

The Purchaser is entitled to free technical support for any ADI product for one year from its date of purchase. Our technical support staff can provide advice concerning installation and operation of ADI products. Services outside of this may incur a charge. Technical support staff will not provide experimental protocols or procedural instructions for conducting experiments. However, information of this type may be provided in the supplied product documentation, or on ADI web sites.

Inquiries

For additional information or service inquiries please contact the nearest ADInstruments office or Authorized Distributor. For contact details see www.ADInstruments.com

Copyright © ADInstruments NZ Ltd, 2022. All rights reserved. PowerLab, LabChart and ADInstruments are registered trademarks of ADInstruments NZ Ltd. Windows 8, Windows 7, Windows Vista and .NET Framework are trademarks of Microsoft Corporation. Apple, the Apple logo, MacOS, and Macintosh are trademarks of Apple Computer Inc. registered in the U.S. and other countries. Acrobat and Adobe are registered trademarks of Adobe Systems Incorporated. Igor is a trademark of Wavemetrics Inc. MATLAB is a registered trademark of The MathWorks Inc. Grass is a trademark of Astro-Med Inc. All other trademarks are the property of their respective owners.