

MLA1052 Pressure Gauge

Accessories

The MLA1052 Pressure Gauge and accessories can be used to calibrate conventional fluid-filled transducers or Millar Mikro-Tip® pressure and pressure-volume catheters.

The kit includes:

3 x Syringes (20 ml)	1 x Sphygmomanometer
3 x Touhy borst adaptors	1 x Silicon tubing (1m)
3 x 4-way luer stop-cocks, with male locking rings	2 x Silicon Tubing (0.5 m)



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1052 Pressure Gauge

For use with:

Conventional Fluid-Filled Transducers

Millar Mikro-Tip Pressure Catheters

Millar Mikro-Tip Pressure-Volume Catheters