

## MLA1213 Needle Electrodes for FE136 (3 pk)

### Accessories

A set of three 1.5 mm shrouded sockets to monopolar needle electrodes (29 gauge) lead wires for use with the FE136 Animal Bio Amp. The lead wires are 122 cm in length and the needle electrodes are made from surgical steel and 12 mm in length. They are shipped non-sterile.



#### Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### Ordering Information:

MLA1213 Needle Electrodes for FE136 (3pk)

For use with:  
FE136 Animal Bio Amp