

## MLA145 Student Disposable Respiratory Kit

### Accessories

The MLA145 Student Disposable Respiratory Kit contains a nose clip, disposable bacterial filter and a mouth piece. The kit is intended for single use only. However the mouth piece is reusable after cold sterilization.



#### Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use. Mouth piece is NOT supplied sterile. Cold sterilize before using.

#### Ordering Information

MLA145 Student Disposable Respiratory Kit

For use with:

Nose Clip

Reusable Mouthpiece

Disposable Filter