

MLA240 Metabolic Accessory Kit

Accessories

The MLA240 Metabolic Accessory Kit contains items for use with the ML240 Metabolic System. It contains a MLT1000L Respiratory Flow Head (1000 L/min), MLA1028 Face Mask with 2 way non rebreathing valve, MLA1081 Flow Head Adapter, MLA1013 35 mm ID Tubing Adapter, MLA1015 Breathing Tube, MLA6024 Desiccant Cartridge, an MLA0343 Drying Tube and Thermistor Temperature Sensor. The items can be purchased separately by quoting the individual part number.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA240 Metabolic Accessory Kit

Includes:

MLT1000L Respiratory Flow Head	MLA1013 35 mm ID Tubing Adapter
MLA1015 Breathing Tube	MLA6024 Desiccant Cartridge
MLA1028 Face Mask	MLA0343 Drying Tube
MLA1081 Flow Head Adapter	Silicon Tube (1mm ID, 30cm length) with Luers x 2
Thermistor Temperature Sensor (incorporates MLT415 Nasal Temperature Probe)	