

## **MLA242 Exercise Breathing Consumables Kit**

### *Accessories*

The MLA242 Exercise Breathing Consumables Kit contains all consumable items for collecting expired air. One consumable kit per PowerLab setup is recommended.

Includes:

- 20 x MLA1081 Flow Head Adapters
- 20 x SP0145 Modified Tubing Adapters
- 20 x Mouthpiece One-Way Valves
- 20 x Douglas Bags (120 litres)
- 20 x Valve End-Caps



#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA242 Exercise Breathing Consumables Kit

For use with:

LabTutor Exercise Physiology Experiments