

MLAYDG Dry Earth Strap

Accessories

The MLAYDG Dry Earth Strap is suitable for use as a ground electrode in biopotential measurements. The strap does not require electrode cream and is suitable for attachment to limbs up to 100 mm (4") in diameter (~315 mm circumference). A snap button on the strap is used to interface with MLA2503 or MLA2505 Lead Wires with snap connectors.



Caution

Read "Statement of Intended Use" on our website or in "Getting Started with PowerLab" before use.

Ordering Information:

MLAYDG Dry Earth Strap


For use with:

MLA2340 3 Lead Shielded Bio Amp Cable

MLA2540 5 Lead Shielded Bio Amp Cable

MLA0313 Lead Wires (3 snap on)

MLA0315 Lead Wires (5 snap on)



NOTE: Lead Wires are connected to the top row of signal pins in Shielded Bio Amp Cables.